

## Protein Kinases

<b>PROTEIN KINASES</b>					
<b>Cat No.</b>	<b>Product Name</b>	<b>Product Description</b>	<b>Size A</b>	<b>Size B</b>	<b>Size C</b>
3HPKA 348	14-3-3B	Recombinant Human Tyr-3/Trp-5 Monooxygenase Activation Protein, Beta	5µg	25µg	1mg
3HPKA 238	14-3-3G	Recombinant Human Protein Kinase C Inhibitor Protein-1	5µg	25µg	1mg
3HPKA 253	14-3-3E	Recombinant Human Tyr-3/Trp-5 Monooxygenase Activation Protein, Epsilon	5µg	25µg	1mg
3HPKA 347	14-3-3 ETA	Recombinant Human Tyr-3/Trp-5 Monooxygenase Activation Protein	2µg	10µg	1mg
3HPKA 254	14-3-3T	Recombinant Human Tyr-3/Trp-5 Monooxygenase Activation Protein, Tau	5µg	25µg	1mg
3HPKA 257	14-3-3Z	Recombinant Human Tyr-3/Trp-5 Monooxygenase Activation protein, Zeta	5µg	20µg	1mg
3HPKA 346	AK3L1	Recombinant Human Adenylate Kinase 3	2µg	10µg	1mg
3HPKA 305	AKT-1	Recombinant Human AKT1	2µg	10µg	1mg
3HPKA 306	ATF-2	Recombinant Human Activating Transcription Factor-2	2µg	10µg	1mg
3HPKA 325	CDK4	Recombinant Human cyclin-dependent kinase 4	2µg	10µg	1mg
3HPKA 310	CDK8	Recombinant Human cyclin-dependent kinase 8	2µg	10µg	1mg
3HPKA 323	C-JUN	Recombinant C-JUN	2µg	10µg	1mg
3HPKA 341	p16-INK4a	Recombinant Human Cyclin-Dependent Kinase Inhibitor 2A	5µg	20µg	1mg
3HPKA 311	CCNC	Recombinant Human cyclin-C	2µg	10µg	1mg
3HPKA 345	CIB1	Recombinant Human Calcium and Integrin Binding 1	2µg	10µg	1mg
3HPKA 307	CREB	Recombinant Human Camp Responsive Element Binding Protein	2µg	10µg	1mg
3HPKA 317	C-SRC	Recombinant Human C-SRC	2µg	10µg	1mg
3HPKA 344	EGFR Sf9	Recombinant Human Epidermal Growth Factor Receptor-Sf9 (ErbB1)	2µg	10µg	100µg
3HPKA 335	EGFR Sf9	Recombinant Human Epidermal Growth Factor Receptor-Sf9 (ErbB1)	1µg	5µg	20µg
3HPKA 336	EGFR	Recombinant Human Epidermal Growth Factor Receptor (ErbB1)	2µg	5µg	10µg
3HPKA 340	c-erbB2	Recombinant Human c-erbB2	2µg	5µg	10µg
3HPKA 343	ErbB2	Recombinant Human Tyrosine Kinase ErbB-2	5µg	20µg	1mg
3HPKA 342	ErbB3	Recombinant Human Tyrosine Kinase ErbB-3	5µg	20µg	1mg
3HPKA 308	EIF2-a	Recombinant Human Eukaryotic Initiation Factor 2-alpha	2µg	10µg	1mg
3HPKA 309	ELK-1	Recombinant Human ELK-1	2µg	10µg	1mg
3HPKA 312	FAK1	Recombinant Human Active Focal Adhesion Kinase	1µg	5µg	10µg
3HPKA 230	FGFR1	Recombinant Human Fibroblast Growth Factor Receptor-1	2µg	10µg	1mg
3HPKA 231	FGFR2	Recombinant Human Fibroblast Growth Factor Receptor-2	2µg	10µg	1mg
3HPKA 232	FGFR3	Recombinant Human Fibroblast Growth Factor Receptor-3	2µg	10µg	1mg
3HPKA 233	FGFR4	Recombinant Human Fibroblast Growth Factor Receptor-4	2µg	10µg	1mg
3HPKA 237	IGF-IR	Recombinant Human Insulin Like Growth Factor-I Receptor	2µg	5µg	10µg
3HPKA 327	bPI3Ka	Recombinant bovine Phosphoinositide 3-kinase alpha p110a/p85a	1µg	2.5µg	5µg
3HPKA 328	bP85a	Recombinant bovine Phosphoinositide 3-kinase a, regulatory subunit	2µg	10µg	20µg

3HPKA 329	PI3K <sup>His</sup>	Recombinant Human Phosphoinositide 3-kinase p110g	1µg	5µg	10µg
3HPKA 330	PI3K <sup>GST</sup>	Recombinant Human Phosphoinositide 3-kinase p110d/p85a	1µg	5µg	10µg
3HPKA 331	PI3K <sup>γD946</sup> <sup>GST</sup>	Recombinant Human Phosphoinositide 3-kinase p110g inactive mutant	1µg	2.5µg	5µg
3HPKA 332	PI3K <sup>His-GST</sup>	Recombinant Human Phosphoinositide 3-kinase p110g <sup>His</sup> /p101 <sup>GST</sup>	1µg	5µg	10µg
3HPKA 333	PI3K <sup>β</sup>	Recombinant Human Phosphoinositide 3-kinase beta P110 <sup>β</sup> /p85a	1µg	2.5µg	5µg
3HPKA 200	PRKACA	Recombinant c-AMP dependant Protein Kinase A catalytic subunit alpha	1µg	5µg	100µg
3HPKA 201	PRKA1RA	Recombinant c-AMP dependant Protein Kinase A regulatory subunit-1 alpha	1µg	3µg	100µg
3HPKA 202	cfPKAR-1a	Recombinant Protein Kinase A regulatory subunit-1 alpha	5µg	25µg	50µg
3HPKA 203	PKACa2/RIa2	Recombinant Protein Kinase A Inactive holoenzyme type 1 alpha	1µg	3µg	100µg
3HPKA 204	PKAr-1a	Recombinant c-AMP dependant Protein Kinase A regulatory subunit 2 alpha	1µg	3µg	100µg
3HPKA 205	PKACa2/RIIa2	Recombinant Protein Kinase A Inactive holoenzyme type 2 alpha	1µg	3µg	100µg
3HPKA 206	PKAkt1/PKBa	Recombinant Human Protein Kinase Akt1/PKB alpha Active Enzyme	5µg	20µg	40µg
3HPKA 207	PKAkt1/PKBa	Recombinant Human Protein Kinase Akt1/PKB alpha Inactive Enzyme	5µg	20µg	40µg
3HPKA 208	CK2a	Recombinant Human Casein Kinase 2 alpha	3µg	10µg	20µg
3HPKA 225	CK2a His	Recombinant Human Protein kinase Casein Kinase 2 alpha, His Tag	2µg	10µg	100µg
3HPKA 209	CK2b Active	Recombinant Human Casein Kinase 2 beta Active	3µg	10µg	20µg
3HPKA 223	CK2b	Recombinant Human Casein Kinase 2 beta	2µg	10µg	100µg
3HPKA 210	ZMCK2a	Recombinant Zea Mays Casein Kinase 2 alpha	3µg	10µg	20µg
3HPKA 211	CK2h	Recombinant Human Casein Kinase 2 Holoenzyme	3µg	10µg	20µg
3HPKA 212	ERK1/MAPK3	Recombinant Human ERK1/MAPK3 Active Enzyme	3µg	10µg	20µg
3HPKA 213	ERK1/MAPK3	Recombinant Human ERK1/MAPK3 Inactive Enzyme	3µg	10µg	20µg
3HPKA 214	ERK2/MAPK1	Recombinant Human ERK2/MAPK1 Active Enzyme	1µg	5µg	20µg
3HPKA 215	ERK2/MAPK1	Recombinant Human ERK2/MAPK1 Inactive Enzyme	3µg	10µg	20µg
3HPKA 226	HK-1	Recombinant Human Hexokinase-1	2µg	10µg	1mg
3HPKA 227	HK-2	Recombinant Human Hexokinase-2	2µg	10µg	1mg
3HPKA 229	HK-3	Recombinant Human Hexokinase-3	2µg	10µg	1mg
3HPKA 236	GCK	Recombinant Human Hexokinase-4	2µg	10µg	1mg
3HPKA 216	JNK2/SAPK1	Recombinant Human JNK2/SAPK1	1µg	5µg	200µg
3HPKA 220	JNK2a-1	Recombinant Human Jun N-terminal kinase 2, alpha 1	2µg	10µg	1mg
3HPKA 314	JNK2a-2	Recombinant Human Jun N-terminal kinase 2, alpha 2	2µg	10µg	1mg
3HPKA 250	LYVE-1 Sf9	Recombinant Human Lymphatic Vessel Endothelial Hyaluronic Acid Receptor 1, Sf9	2µg	10µg	1mg
3HPKA 252	LYVE-1	Recombinant Human Lymphatic Endothelial Hyaluronan Receptor-1	2µg	5µg	10µg
3HPKA 251	mLYVE-1 Sf9	Recombinant Mouse Lymphatic Vessel Endothelial Hyaluronic Acid Receptor 1, Sf9	2µg	10µg	1mg
3HPKA 256	MEK1	Recombinant Human Mitogen Activated Kinase 1	1µg	5µg	10µg
3HPKA 339	PKM2	Recombinant Human Tumour Type M2 Pyruvate Kinase	2µg	10µg	1mg
3HPKA 337	NCK2	Recombinant Human NCK2	2µg	5µg	10µg
3HPKA 217	p38a/SAPK2a	Recombinant Human p38	1µg	5µg	20µg

		alpha/SAPK2 alpha			
3HPKA 218	PKC-a	Recombinant Human Protein Kinase C alpha	2µg	10µg	1mg
3HPKA 318	PKC-b1	Recombinant Human Protein Kinase C beta 1	2µg	10µg	1mg
3HPKA 319	PKC-b2	Recombinant Human Protein Kinase C beta 2	2µg	10µg	1mg
3HPKA 320	PKC-g	Recombinant Human Protein Kinase C gamma	2µg	10µg	1mg
3HPKA 321	PKC-d	Recombinant Human Protein Kinase C delta	2µg	10µg	1mg
3HPKA 322	PKC-e	Recombinant Human Protein Kinase C e	2µg	10µg	1mg
3HPKA 324	PKC-q	Recombinant Human Protein Kinase C q	2µg	10µg	1mg
3HPKA 228	PPAR-g	Recombinant Human Peroxisome Proliferator Activated Receptor Gamma	2µg	10µg	1mg
3HPKA 334	PPAR-g (1-477)	Recombinant Human Peroxisome Proliferator Activated Receptor Gamma (1-477)	2µg	5µg	10µg
3HPKA 224	PSPH	Recombinant Human Phosphoserine Phosphatase	5µg	20µg	1mg
3HPKA 219	PTP-1B	Recombinant Human Protein Tyrosine Phosphatase 1B	10µg	50µg	1mg
3HPKA 222	PP2Ca	Recombinant Human Protein Phosphatase 2C alpha	10µg	50µg	1mg
3HPKA 221	SHP-1	Recombinant Human SHP-1	5µg	20µg	1mg
3HPKA 315	STAT1	Recombinant Human STAT1	2µg	10µg	1mg
3HPKA 316	STAT3	Recombinant Human STAT3	2µg	5µg	10µg
3HPKA 246	TIE-1 Fc	Recombinant Human TIE-1 Fc Chimera	2µg	10µg	1mg
3HPKA 247	TIE-1 Fc	Recombinant Mouse TIE-1 Fc Chimera	2µg	10µg	1mg
3HPKA 248	TIE-2 Fc	Recombinant Human TIE-2 Fc Chimera	2µg	10µg	1mg
3HPKA 249	TIE-2 Fc	Recombinant Mouse TIE-2 Fc Chimera	2µg	10µg	1mg
3HPKA 255	TPMT	Recombinant Human Thiopurine S-methyltransferase	5µg	20µg	1mg
3HPKA 240	VEGFR1	Recombinant Human Vascular Endothelial Growth Factor receptor-1	2µg	10µg	100µg
3HPKA 338	VEGFR1 D3 His	Recombinant Human Vascular Endothelial Growth Factor receptor-1 D1-3, His Tag	2µg	10µg	1mg
3HPKA 234	VEGFR1 D3	Recombinant Human Vascular Endothelial Growth Factor receptor-1 D1-3	2µg	10µg	100µg
3HPKA 235	VEGFR1 D4	Recombinant Human Vascular Endothelial Growth Factor receptor-1 D1-4	2µg	10µg	100µg
3HPKA 239	VEGFR1 D5	Recombinant Human Vascular Endothelial Growth Factor receptor-1 D1-5	2µg	10µg	100µg
3HPKA 241	VEGFR1 D7	Recombinant Human Vascular Endothelial Growth Factor receptor-1 D1-7	2µg	10µg	1mg
3HPKA 242	VEGFR2	Recombinant Human Vascular Endothelial Growth Factor receptor-2	2µg	10µg	1mg
3HPKA 243	VEGFR2 Fc	Recombinant Human Vascular Endothelial Growth Factor receptor-2 Fc Chimera	2µg	10µg	1mg
3HPKA 244	VEGFR3	Recombinant Human Vascular Endothelial Growth Factor receptor-3	2µg	10µg	1mg
3HPKA 245	VEGFR3 Fc	Recombinant Human Vascular Endothelial Growth Factor receptor-3 Fc Chimera	2µg	10µg	1mg
3HCKI 268	CKBB	Recombinant Human Creatine Kinase BB Isoenzyme	10µg	50µg	1mg
3HCKI 269	CKMB Type-1	Recombinant Human Creatine Kinase MB Isoenzyme Type-1	10µg	50µg	1mg
3HCKI 270	CKMB Type-2	Recombinant Human Creatine Kinase MB Isoenzyme Type-2	10µg	50µg	1mg
3HCKI 273	CKMM	Human Creatine Kinase MM	200µg	1mg	10mg
3HCKI 271	CKMM Type-1	Recombinant Human Creatine Kinase MM Isoenzyme Type-1	10µg	50µg	1mg



3H Biomedical AB  
Dag Hammarskjölds väg 34A  
SE-751 83 Uppsala, Sweden  
Phone: +46 (0)18 50 44 00  
Fax: +46 (0)18 51 38 44  
E-mail: info@3hbiomedical.com

3HCKI 272	CKMM Type-3	Recombinant Human Creatine Kinase MM Isoenzyme Type-3	10µg	50µg	1mg
-----------	-------------	--	------	------	-----